

NEW PX SLIM™

The Penumbra Coil 400™ Embolization System

Exceptionally Stable Occlusion



Frame/Fill

Fill

Fill/Finish



Complex Standard

(3–32 mm)



Complex Soft



J Soft

(up to 25 cm)



Complex Extra Soft

(2–4 mm)





Curve Extra Soft

Penumbra 

The Penumbra Coil 400™ Embolization System

Frame	Fill	Fill/Finish	
COMPLEX STANDARD	COMPLEX SOFT	COMPLEX EXTRA SOFT	CURVE EXTRA SOFT



Catalog Number	Secondary Diameter (mm)	Length (cm)
4002C0306	3	6
4002C0310	3	10
4002C0408	4	8
4002C0412	4	12
4002C0506	5	6
4002C0510	5	10
4002C0610	6	10
4002C0615	6	15
4002C0710	7	10
4002C0715	7	15
4002C0720	7	20
4002C0815	8	15
4002C0820	8	20
4002C0830	8	30
4002C0915	9	15
4002C0925	9	25
4002C0935	9	35
4002C1020	10	20
4002C1030	10	30
4002C1040	10	40
4002C1135	11	35
4002C1145	11	45
4002C1235	12	35
4002C1245	12	45
4002C1335	13	35
4002C1348	13	48
4002C1450	14	50
4002C1557	15	57
4002C1660	16	60
4002C1857	18	57
4002C2060	20	60
4002C2260	22	60
4002C2457	24	57
4002C2860	28	60
4002C3260	32	60


Catalog Number	Secondary Diameter (mm)	Length (cm)
4004C0303	3	3
4004C0304	3	4
4004C0306	3	6
4004C0310	3	10
4004C0404	4	4
4004C0406	4	6
4004C0408	4	8
4004C0412	4	12
4004C0509	5	9
4004C0513	5	13
4004C0610	6	10
4004C0615	6	15
4004C0620	6	20
4004C0710	7	10
4004C0715	7	15
4004C0720	7	20
4004C0815	8	15
4004C0820	8	20
4004C0830	8	30
4004C0925	9	25
4004C0935	9	35
4004C1020	10	20
4004C1030	10	30
4004C1040	10	40
4004C1135	11	35
4004C1145	11	45
4004C1235	12	35
4004C1245	12	45
4004C1335	13	35
4004C1348	13	48
4004C1440	14	40
4004C1545	15	45
4004C1645	16	45

Catalog Number	Secondary Diameter (mm)	Length (cm)
4006C0202	2	2
4006C0203	2	3
4006C0204	2	4
4006C0303	3	3
4006C0304	3	4
4006C0305	3	5
4006C0306	3	6
4006C0308	3	8
4006C0310	3	10
4006C0404	4	4
4006C0406	4	6
4006C0408	4	8
4006C0410	4	10

Catalog Number	Secondary Diameter (mm)	Length (cm)
4006U0201	2	1
4006U0202	2	2
4006U0203	2	3
4006U0204	2	4
4006U0302	3	2
4006U0303	3	3
4006U0304	3	4
4006U0305	3	5
4006U0306	3	6
4006U0308	3	8
4006U0310	3	10
4006U0404	4	4
4006U0406	4	6
4006U0408	4	8
4006U0410	4	10


Fill

J SOFT




Catalog Number	Length (cm)
4004J7	7
4004J10	10
4004J15	15
4004J25	25

Penumbra Coil Detachment Handle



Catalog Number	Description
DH1	Detachment Handle
DH5	Detachment Handle (5 pack)

PX SLIM™ Delivery Microcatheter



Catalog Number	Tip Shape	ID	Distal and Proximal OD	Guidewire Compatibility	Tip Markers	Length
PXS LIMSTR	Straight	0.025"	2.6/2.95F	≤0.020"	2	150 cm
PXS LIM045	45°	0.025"	2.6/2.95F	≤0.020"	2	150 cm
PXS LIM090	90°	0.025"	2.6/2.95F	≤0.020"	2	150 cm
PXS LIM130	130°	0.025"	2.6/2.95F	≤0.020"	2	150 cm



www.penumbrainc.com

Penumbra, Inc. USA
1351 Harbor Bay Pkwy
Alameda, CA 94502
USA
1.888.BRAIN.06
1.888.272.4606
T 1.510.748.3200
F 1.510.748.3232
order@penumbrainc.com
info@penumbrainc.com

Penumbra Europe GmbH
Am Borsigturm 44
13507 Berlin
Germany
T +49 30 2005 676-0
F +49 30 2005 676-10
order@penumbrainc.de
info@penumbrainc.de

Penumbra Neuro Australia PTY LTD
Level 10, 1 Margaret Street
Sydney NSW 2000
Australia
T +61-1300 817 025
F +61-1300 817 026
order.anz@penumbrainc.com